

Mammomat Select is not commercially available in all countries. Due to regulatory reasons, its future availability cannot be guaranteed. Please contact your local Siemens organization for further information.

On account of certain regional limitations of sales rights and service availability, we cannot guarantee that all products/ services/features included in this brochure are available through the Siemens sales organization worldwide. Availability and packaging may vary by country and are subject to change without prior notice.

The information in this document contains general descriptions of the technical options available and may not always apply in individual cases.

The product names and/or brands referred to are the property of their respective trademark holders.

Siemens reserves the right to modify the design and specifications contained herein without prior notice. Please contact your local Siemens sales representative for the most current information.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources and waste conservation), we may recycle certain components where legally permissible. For recycled components we use the same extensive quality assurance measures as for factory-new components.

Any technical data contained in this document may vary within defined tolerances. Original images always lose a certain amount of detail when reproduced.

SIEMENS



Mammomat Select


Makes your day the easy way

[siemens.com/mammomat-select](https://www.siemens.com/mammomat-select)

Siemens Healthcare Headquarters
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
[siemens.com/healthcare](https://www.siemens.com/healthcare)

Order No. A91XP-30008-1C2-7600 | Printed in Germany | CC XP 3440 02162. | © Siemens Healthcare GmbH, 2016

[siemens.com/healthcare](https://www.siemens.com/healthcare)

A Siemens Mammomat Select mammography machine is shown from a low-angle perspective. The machine is white with a pink control panel on top. The Siemens logo is visible on the upper part of the machine. The text "Mammomat Select" is printed on the side. The machine is positioned in a room with large windows in the background.

Mammomat Select Makes your day the easy way

Providing access to mammography and therefore to early breast cancer detection can be a challenge when budgets are limited. An analog system that is easy to use and provides images with the right dose, is a cost-effective way to provide women with professional care.

Mammomat Select* is a system that is the outcome of more than 30 years of expertise in analog mammography. Enhanced with smart features for ease of use, the right dose and low total costs of ownership, it will help you give more women access to high-quality mammography .

In short, with Mammomat Select, you are able to offer more women worldwide access to mammography in keeping the required standard of care.

Easy to touch screen

Features like EasyTouch and isocentric rotation allow you to streamline your workflow.

Easy to dose right

Personalized OpDose with Automatic Exposure Control makes individually selecting the right dose for each woman very simple.

Easy to invest low

Mammomat Select comes with an affordable price tag, low service costs and enables high patient throughput.

Discover how these attractive features can significantly improve your daily screening and diagnostics. Because detail matters in mammography.

*Mammomat Select is not commercially available in all countries.
Due to regulatory reasons its future availability cannot be guaranteed.
Please contact your local Siemens organization for further information.

Easy to touch screen

Reduce staff training requirements to a minimum with the EasyTouch user interface.

Clear and compact controls

The intuitive EasyTouch allows full access to the exposure controls from a single compact console, either with the integrated radiation shield or from a separate control room.

Automated procedures

EasyTouch comes with intuitive symbols and the next workflow step is highlighted, allowing for easy operation throughout the entire mammography exam. Individual changes for each patient, such as chamber selection, are possible at the touch of a button.

User-defined configuration tabs

You have full access to the system configuration and can directly change the dose settings, exposure light levels, compression force and more with just one touch.

Easy to dose right

Give each and every patient the same high standard of care with dose saving features that keep exposure to a minimum.

OpDose with one touch

OpDose can be selected on the intuitive EasyTouch with a single tap. It works quickly and easily in every position and view.

The right dose for every woman

OpDose uses breast thickness and Automatic Exposure Control (AEC) measurements to achieve optimal image quality for each woman at the right dose.

*Option

Easy to invest low

Reduce the total cost of ownership with a system with compact space requirements, easy installation and flexible service packages.

High throughput

Isocentric rotation and fully motorized system movements enable fast positioning, and, combined with the efficient EasyTouch, let you increase patient throughput.

Flexible configuration

The compact design and integrated generator reduce space requirements. The ability to configure the system with a Mo/Mo or Mo/Rh anode-filter combination* and an 18 x 24 or 24 x 30 cm bucky table allows you to install the system according to your needs.

Expanded applications

Mammomat Select offers options such as magnification, spot imaging, 2D biopsy and a range of compression plates for flexible positioning. The interfaces for a Film ID camera or CR reader further increase the scope of possible applications.

*Option



EasyTouch allows easy operation at the touch of a button.

Join the revolution: because detail matters.

Mammomat Select is a system that provides access to mammography and therefore to early breast cancer detection, which can be a challenge when budgets are limited.

It combines proven expertise in analog mammography, enhanced with smart features for ease of use, low dose and cost of ownership to improve your standard of care.

Features like EasyTouch and Isocentric rotation allow you to streamline your workflow. Personalized OpDose with Automatic Exposure Control makes individually selecting the right dose for each woman very simple. And the system's low operating costs, attractive service packages and high patient throughput also help you focus on financial performance.

Mammomat Select is easy to use and provides images with the right dose. It is a cost-effective way in providing women the standard of care they require.

Service tailored for you

Siemens offers a range of comprehensive service options that can be precisely tailored to your requirements, helping you to optimize system utilization and improve your workflow. Ranging from preventive maintenance to technical phone support and system updates, our Performance and Protect Plans are proven service packages that include proactive services and rapid technical on-site responses to help improve reliability while offering predictable costs, low risk and great efficiency.



Facts

How much mammography every day?

3,000,000,000

women live on the planet.

650,000,000

are eligible for breast cancer screening programs.

225,000,000

mammography exams per year.

615,000

women need easy access to mammography every day.

Makes your day the easy way.



Fully motorized system movements
enable quick positioning – at the touch of a button

Personalized OpDose with AEC
individually selects the right dose for every woman

Seamless isocentric rotation
enables rapid positioning for the next view without changing detector height

Compact installation
and reduced space requirements due to integrated generator

Makes your day the easy way

Tube with anode-filter combination
lets you choose between a Mo/Mo or Mo/Rh combination*

Flexible configuration
with an 18 x 24 cm or 24 x 30 cm bucky table as well as magnification, spot and 2D biopsy options

EasyTouch
allows easy operation at the touch of a button

*Option